LENIRE®

Neuromod Devices Ltd., Rainsford St., Dublin, D08 R2YP, Ireland.
www.neuromoddevices.com

Lenire Tinnitus Treatment System. User Manual

- DE Lenire-Tinnitus-Behandlungssystem. Benutzerhandbuch
- FR Manuel de l'utilisateur Système de traitement des acouphènes
- IT Lenire Manuale d'uso Lenire Tinnitus Treatment System
- NL Lenire Tinnitus Behandelingssysteem. Gebruikershandleiding

IFU-0026 V1.0

English *page 2*

Using this manual

This manual is supplied with the Lenire system to provide you with information on the product, treatment regime and instructions for its safe and effective use. You should read and understand these instructions before use.

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Table of Contents

Using this manual	1
Section A - General information	3
Chapter 1 - Introduction	3
Chapter 2 - Intended use	3
Chapter 3 - Contraindications	3
Chapter 4 - Risks associated with use	3
Chapter 5 - General warnings	4
Chapter 6 - Product and package labelling	5
Section B - Using Lenire	7
Chapter 7 - Product components and accessories	7
Chapter 8 - Controls and indicators	8
Chapter 9 - Controller port	11
Chapter 10 - Charging the Controller and Headphones	11
Chapter 11 - Device setup	13
Chapter 12 - Tonguetip setup	14
Chapter 13 - Starting and controlling a treatment session	15
Chapter 14 - Treatment guidelines	17
Section C - Care, storage maintenance and disposal	18
Chapter 15 - Cleaning	18
Chapter 16 - Storage and transport	19
Chapter 17 - Maintenance	20
Chapter 18 - Disposal	20
Section D - Troubleshooting	21
Chapter 19 - Frequently asked questions	21
Chapter 20 - Bluetooth pairing / un-pairing	23
Chapter 21 - Servicing	24
Chapter 22 - Device errors	24
Section E - Technical information	25
Chapter 23 - Environmental operating conditions	25
Chapter 24 - Radio frequency (RF) specifications	25
Chapter 25 - Guidance and declaration for electromagnetic compatibility	25
Chapter 26 - Materials in direct contact with human tissue	26
Chapter 27 - Declaration of conformity radio equipment directive 2014/53/EU	26
Section F - Customer support contact details	27

Section A

General information

This section outlines the general information relating to the Lenire® system, including important safety information.

Chapter 1 - Introduction

Lenire is a medical device intended to reduce the symptoms of tinnitus. It comprises a handheld Controller and a Tonguetip® which, delivers gentle, electrical stimulation to the tongue, and is used in conjunction with wireless Headphones that deliver audio stimulation. A Charger is included that can be used to charge both the Controller and the Headphones.

The audio and tongue stimulation are configured and calibrated to your individual characteristics during the initial fitting procedure by a trained clinician.

Chapter 2 - Intended use

The Lenire device is intended to be used by tinnitus sufferers of at least 18 years of age to alleviate the symptoms of chronic, subjective tinnitus. It should be used for 30 to 60 minutes per day, every day, for at least 10 weeks, after which you may experience sustained alleviation of your tinnitus symptoms.

The Lenire device can be used in the home environment after consultation and fitting by a suitably qualified clinician.

Lenire is intended for prescription use only. Use by a person for whom the device has not been configured may cause discomfort, or temporary exacerbation of tinnitus loudness.

Chapter 3 - Contraindications

Caution! – Do not use the device if:

- You have a pacemaker, defibrillator or any other active implantable device, unless directed by a doctor.
- You are pregnant, unless directed by a doctor.
- You suffer from epilepsy or any other condition that may result in loss of consciousness.
- You suffer from any condition that causes impaired sensitivity of the tongue.
- You have lesions, sores, or inflammation of the oral cavity, unless directed by a doctor.

Chapter 4 - Risks associated with use

Potential risks associated with use of the Lenire device include

Contact between the Tonguetip and dental retainers or metal fillings may cause temporary discomfort.

There is a possibility that some people may experience a mild allergic reaction to the materials in the device.

You may experience a mild numbing sensation or a discomfort on the tongue during treatment. If this occurs, reduce the tongue stimulation. If the symptoms persist, contact your clinician.

Use of the treatment may result in mild fluctuations in tinnitus symptoms, such as tinnitus loudness.

Headaches may accompany use of the treatment device.

Chapter 5 - General warnings

Do not use the Tonguetip when going to sleep or under the influence of alcohol or other substances that cause drowsiness.

Keep the Tonguetip and Charger out of the reach of children, to avoid possible choking and strangulation hazards.

To avoid the risk of cross-contamination, the Tonguetip should not be shared with any other person, as it is a single-patient, multi-use device.

Do not use the Tonguetip if you have lesions, sores, or inflammation of the oral cavity. Discontinue use of the device until symptoms have resolved.

Do not leave the device in the sun, near an open fire or on a radiator.

The use of rechargeable lithium ion batteries is very common and they are generally considered safe to use. With any energy storage device, they carry safety risks including overheating. Only use the charger provided with your device and always follow the instructions in this manual when charging your Lenire device.

Do not attempt to replace the Lithium Ion battery in your Lenire device. Replacement of the battery should be performed by a qualified individual only; attempts to replace the battery by an inadequately trained person may result in risks including overheating and fire.

Do not use in conjunction with a hearing aid.

To avoid the risk of a burn do not puncture, crack or crush the Controller.

Chapter 6 - Product and package labelling

The label on the Controller, the Tonguetip plug, Carry Case and the product packaging provide important information. The meanings of the associated symbols are described below, and in the context of the relevant sections throughout this user manual.

Use-by date (year and month). This is the date after which the device should no longer be used.

Do not use if the packaging is damaged or previously opened.

Consult the user manual.

This symbol indicates the need for the user to consult the User Manual for warnings and precautions associated with this medical device.

Manufacturer's details (name and address) and date of manufacture (year and month).

The device's serial number, model number and lot number appears next to these symbols. You will need these numbers if you call customer support for assistance. See page 53 for the manufacturer's contact details.

IP22

Ingress Protection (IP) Rating: IP ratings describe how effectively an enclosure is sealed.

IP22 rated equipment is protected from objects greater than 12mm in diameter being inserted e.g. fingers, and against small amounts of vertically dripping water.

Non-ionizing radiation.

This symbol appears next to the model number of the Charger which should be used to charge the Controller.

This symbol signifies that the device should not be disposed of as domestic waste. This ensures compliance with the Battery Directive 2006/66/EC.

The logo for this product range (Lenire tinnitus treatment system).

This symbol (direct current) identifies a charging port.

This symbol represents the Tonguetip.

The USB (Universal Serial Bus) symbol – this is the connection that is used by the Clinician when they configure your device.

This symbol indicates that the Tonguetip provides a high degree of protection against electric shock (Type BF Applied Part).

Storage and transport temperature limits. The lower limit is indicated near the lower horizontal line and the upper limit is indicated near the upper line.

Storage and transport humidity limits. The lower limit is indicated near the lower horizontal line and the upper limit is indicated near the upper line.

Keep dry.

Fragile, handle with care.

This symbol, the CE Mark, signifies compliance with applicable EU Directives.

The number (2797) is the Notified Body number for the British Standards Institute (BSI), who approved Lenire for CE Marking.



Section B

Using Lenire

This section describes the components, controls and indicators of Lenire and outlines the steps required to prepare for and commence treatment.

Chapter 7 - Product components and accessories

In the Lenire system box, you will find:

Controller	The Lenire Controller (Model Number CR-201) is used to “control” the treatment session i.e. you can start, pause and resume the treatment session, adjust the volume of the sound, adjust the level of the tongue stimulation, and turn the device on or off.	
Tonguetip	The Tonguetip (Model Number TP-101) is a device which provides gentle, electrical stimulation to the tongue.	
Charger	The Charger (XP Power model VEP08US05-XE0767) is used to charge the Controller and the Headphones.	
Carry Case	The Carry Case (Model Number CE-401) is a portable case which can be used to store the Controller and Tonguetip in-between uses.	

The Lenire system box also contains a Quick Start Guide which outlines the basic steps needed to commence a treatment session.

7.1 - Headphones

Lenire must be used in conjunction with a set of Bluetooth Headphones. For best performance, it is essential that you use AKG K845BT Headphones.

These Headphones are an accessory that must be used with the Lenire system and are made available to you at the time of purchase.

Lenire is designed to work only with the AKG K845 BT Headphones, and will not connect to any other Bluetooth devices.

7.2 - Batteries

The Controller is powered by a 3.7V Lithium Ion battery. The battery is charged using the supplied charger. See chapter 10 for instructions on how to charge the Controller. The Controller battery is pre-fitted and is not intended to be replaced.

The Headphones are powered by 3.7V Lithium Ion battery. The Headphones are also charged using the supplied charger.


See chapter 10 for instructions on charging the Controller and Headphones.

If there are difficulties charging the Controller or the Headphones, see the section on Troubleshooting or contact customer support (see section F).

Chapter 8 - Controls and indicators

8.1 - Treatment controls

There are three control buttons on the front of the Controller:

1. Power button

Hold this button for 2 seconds to switch the device on.

Hold this button for 2 seconds to switch the device off.

2. Play button


Press this button to start or resume the treatment session.

3. Pause button

Press this button to pause the treatment session.

8.2 - Indicators

The status of the device is indicated by several lights surrounding the buttons on the top of the Controller.

Battery indicator

Green	The battery power is ok.
Orange	The battery has less than 1 hour of power remaining. Please charge the device before next use (see page 11).
Red	The battery is exhausted. The device will shut off automatically. Please charge the device before next use (see page 11).
Orange (pulsing)	The Charger is connected and the battery is charging.
Green (pulsing)	The Charger is connected and the battery is fully charged.

Headphones indicator

Blue (pulsing)	The Controller is trying to establish a wireless connection to the Headphones.
Green	The Controller is connected wirelessly to the Headphones.
Blue (pulsing rapidly)	The Controller is attempting to pair with the Headphones, see page 23 for details.

Tonguetip indicator

Blue (pulsing)	There is no Tonguetip connected.
Green (pulsing)	The Tonguetip is connected but it is not in contact with the tongue.
Green	The Tonguetip is connected and it is in contact with the tongue.
Orange or orange (pulsing)	The Tonguetip has less than 10 hours of use remaining before expiry. Please purchase a new Tonguetip to avoid disruption to your treatment regime.
Red	The Tonguetip has expired. Please purchase a new Tonguetip to continue treatment.

Treatment indicator

White (pulsing)	The treatment session is paused.
White (animated)	The treatment session is active – the white light rotates around the buttons.
Green (pulsing)	The treatment session has completed.


When the device is powered off, none of the lights will be illuminated.


8.3 - Headphones volume and Tonguetip stimulus controls

Tonguetip stimulation control


The Tonguetip stimulation intensity can be adjusted to a comfortable level by pressing the buttons indicated below. Please see page 16 for guidelines on setting the stimulation intensity.


Indication of default (optimum) Tonguetip stimulus intensity

Indication of maximum Tonguetip stimulus intensity


Indication of minimum Tonguetip stimulus intensity


Headphones volume control


The Headphones volume can be adjusted to a comfortable level by pressing the buttons indicated below. Please see page 17 for guidelines on setting the volume level.


Indication of default (optimum) Headphones volume level

Indication of maximum Headphones volume level

Indication of minimum Headphones volume level


The Headphones volume and Tonguetip stimulus are calibrated to the optimum levels during device fitting by the clinician. For further details, please see page 16.

Chapter 9 - Controller port

The Controller port is used to connect the Tonguetip and the Charger to the Controller. For safety, the Tonguetip and the Charger cannot be connected to the Controller at the same time.

The following symbols are adjacent to the Controller port:

This symbol (direct current) identifies a charging port.

This symbol represents the Tonguetip.

The USB symbol – this is the connection that is used by the clinician when they configure your device.

This symbol indicates that the Tonguetip provides a high degree of protection against electric shock.

The Tonguetip and the supplied Charger are the only equipment that should ever be plugged into the Controller port.

Chapter 10 - Charging the Controller and Headphones

Before Lenire is used for the first time, it is recommended that the Controller and the Headphones be fully charged.

Only use the supplied Charger with the Lenire device.

To avoid the risk of electric shock, do not connect or disconnect the Charger from the mains if it, or your hands, are wet.


The supplied Charger can be used for charging both the Headphones and the Controller



10.1 - Charging the Controller

To charge the Controller:

- Plug the cable end of the supplied Charger (XP Power model VEP08US05-XE0767) into the Controller port of the Controller.
- Plug the Charger into a mains electrical socket.

While charging, the Battery Indicator will pulse orange. Once the battery is fully charged, the Battery Indicator will pulse green.

An empty battery should take no more than 3 hours to fully recharge, and a full charge should allow the Lenire device to be used for at least 8 hours.


Once the battery is fully charged, the Charger should be disconnected to help conserve the life of the battery in the Controller.

10.2 - Charging the Headphones


To charge the Headphones:

- Plug the cable end of the supplied Charger (XP Power model VEP08US05-XE0767) into the USB port of the Headphones.
- Plug the Charger into a mains electrical socket.


While charging, the light on the Headphones will illuminate red. Once the battery is fully charged, the light on the Headphones will extinguish.

Battery charging

Battery charged

When the Headphones battery is running low, the light beside the Headphones power button will blink red and the Headphones will emit a sequence of low-pitch beeps.

An empty battery should take no more than 3 hours to fully recharge, and a full charge should allow the Headphones to be used for at least 8 hours.

Once the battery is fully charged, the Charger should be disconnected to help conserve the life of the battery in the Headphones.

Chapter 11 - Device setup

To commence a treatment session, the Lenire device should be set up as follows:

11.1 - Controller and Headphones setup

The Headphones should be worn so they cover both ears and should be adjusted to maximise comfort before commencing the treatment session.

It is essential that the left and right earphones cover the left and right ears, respectively. This is because the device calibration is unique to each ear – it is based on your audiogram. Failure to put on the Headphones correctly could result in excessive loudness in one ear, which could cause discomfort, or incorrect loudness that may reduce the effectiveness of the treatment.


The **left** and **right** sides of the earphones are marked **L** and **R** respectively on the cloth inside the foam cushions.

11.2 - Connecting the Headphones with the Controller

The Controller connects to the Headphones wirelessly via Bluetooth by following these steps:

- Turn the Headphones on by briefly pressing the power button (the small circular button). About one second later, the blue light beside this button should start to flash blue, or blue and red if the headphones battery is low.
- Turn on the Controller, by pressing the power button for two seconds. The Controller will attempt to connect to the Headphones during which time the Headphones Indicator will pulse blue.
- Once the connection to the Headphones has been established, the Headphones Indicator on the Controller will illuminate green.

The above instruction assumes that the Headphones are already 'paired' with the Controller. This should be the case in most situations, as the clinician will have paired the Headphones during the device fitting. If they are not paired, or if you are having difficulty connecting to the Headphones, then you should follow the instructions on page 23.

During treatment, the Controller can be kept adjacent to you, e.g. in a pocket or a bag, or in another dry location that is close by, e.g. on a bedside table.

Chapter 12 - Tonguetip setup

The Tonguetip is supplied within a protective pouch (made of clear film, designed to be peeled open in the area closest to the V shaped seal).


If the packaging is damaged, the Tonguetip should not be used and it should be returned to the point of sale.

Before connecting and using the Tonguetip for the first time, it should be cleaned as per the instructions on page 18.

The Tonguetip must be connected to the Controller by plugging it into the Controller port, such that the label side of the plug faces away when viewed from the front of the Controller.

Once connected, the Tonguetip Indicator should change from pulsing blue to pulsing green.

The Tonguetip should be placed in your mouth, just behind your upper front teeth (on the top of the tongue towards the front of the mouth) with the electrodes facing down.

As you close your mouth, the electrodes should make contact with the tip of your tongue and rest there comfortably.

Once adequate contact is made with the tongue, the Tonguetip Indicator should stop pulsing (it should illuminate green without pulsing).

For maximum comfort, the mouth and jaw should remain closed but relaxed for the duration of the treatment session.

To avoid the risk of electrical shock or burns, do not connect the controller or Tonguetip to any equipment other than those specified in this user manual.


The Tonguetip model TP-101 is designed to function for 360 treatment sessions (180 hours) before expiring. Indications of imminent expiry and actual expiry are outlined on page 9.

Chapter 13 - Starting and controlling a treatment session

13.1 - Starting a treatment session

As soon as the Headphones and the Tonguetip are connected, which will be indicated by the Headphones Indicator and the Tonguetip Indicator illuminating green, you are ready to start a treatment session.

Once you are in a comfortable position, press the play button to commence the treatment. The Treatment Indicator (white) will rotate around the display.

A single treatment session lasts approximately 30 minutes from the time the play button is pressed.

Towards the end of the treatment session, the sound will gradually fade out over a period of 2 minutes to a lower volume level.

It is essential that you continue to listen to the sound during the fade out period, i.e. until the treatment session is complete.

Once the treatment session is complete, the Treatment indicator will pulse green.

The Controller will power off automatically after 20 minutes of being idle (no treatment session in progress). The Headphones will automatically power off 30 minutes thereafter. It is not necessary to explicitly power off the Controller or the Headphones after the treatment session has completed.

13.2 - Pausing and resuming a treatment session

It is recommended to complete each treatment in a single, continuous 30-minute session to maximise the effectiveness of the treatment.

However, the treatment session can be paused at any time by pressing the Pause Button. The Treatment Indicator will pulse white instead of rotating around the display.

The treatment session will automatically pause if the Tonguetip loses contact with the tongue for more than 30 seconds.

The duration of treatment during a session will remain at 30 active minutes of stimulation. For example, if the treatment is paused for 5 minutes, the total duration of the session would be 35 minutes.

To resume treatment, press the play button while also ensuring that the Tonguetip is in contact with the tongue.

If you do pause a treatment session, it is recommended that you resume again within 5 minutes. It is also recommended that you do not pause the treatment more than once within a single treatment session in order to maximise the effectiveness of the treatment.

13.3 - Adjusting the stimulus levels

13.3.1 - Tonguetip stimulus level control

It is recommended to leave the tongue stimulus level set to the default level unless either:

- There is a strong sensation from the Tonguetip that causes discomfort. In this case the tongue stimulus level can be adjusted downwards to a comfortable level.
- There is a weaker tongue sensation than was perceived during the initial training you received as part of the device fitting. In this case the tongue stimulus level may be adjusted upwards.

The tongue stimulus is calibrated to your sensation thresholds during device fitting by the clinician. The treatment stimulus is then set to be sufficient for optimum effectiveness and may not be perceptible. If the stimulus level is set to the default level (i.e. the level that the device defaults to at the start of each treatment session), then the stimulus is functioning properly as long as the Tonguetip Indicator is illuminated green (i.e. when it is not pulsing).

13.3.2 - Headphones volume control

The Headphones volume level should be left at the default level unless either:

- The treatment environment is noisy, in which case the volume level can be increased.
- The volume is so high that it causes any discomfort, or prevents you from relaxing, in which case the volume level can be decreased.

The Headphones volume levels are calibrated to your hearing thresholds for optimum effectiveness.

Depending on your audiological profile, you may perceive the sound levels to be soft. It is important that you seek a quiet environment for your treatment sessions.

The volume control buttons on the Headphones should not be used to adjust the volume. Please only use the volume control buttons on the Controller.

Depending on your hearing profile, the range over which you can increase the Headphones volume level may be limited. This is to ensure that the volume levels remain within safe limits, so as not to cause hearing damage.

Chapter 14 - Treatment guidelines

The following guidelines should be followed to maximise the effectiveness of the treatment programme:

- While the treatment is intended for use between 30 to 60 minutes per day, it is recommended that you try to complete at least a full 30 minutes at each daily treatment.
- Try to keep the tongue stimulation level and the Headphones volume level at the default settings, as these are the optimum levels calibrated to your individual characteristics.
- Try to complete each treatment session in a place that is quiet, relaxing and comfortable. Focusing your attention away from distractions such as the TV may improve the effectiveness of the treatment.
- Try to keep a routine for the treatment sessions, where the treatment is delivered at the same time(s) each day.
- It is not unusual for tinnitus symptoms to disimprove or fluctuate. It is recommended that you continue treatment even during periods where the symptoms disimprove, for the 10 weeks minimum recommended treatment duration. However, if there is a significant disimprovement in your symptoms it is recommended to seek advice from your clinician.

Do not use whilst showering, in rain, or in a moist atmosphere such as a steam bath or sauna. Do not use whilst applying items such as, perfume, cosmetics, after-shave, hairspray or suntan lotion.

Section C

Care, storage, maintenance and disposal

This section details the correct procedures for cleaning, storage, maintenance and disposal of the Lenire system.

Chapter 15 - Cleaning

It is advisable to clean the Controller, the Headphones, the Carry Case and the Tonguetip after prolonged disuse or storage, or after exposure to conditions that could affect their hygiene.

15.1 - Instructions for cleaning of the Controller, the Headphones and the Carry Case

Before cleaning the Controller or Headphones, ensure that they are switched off and not connected to the Charger.

- The exterior surfaces of the Controller, the Headphones and Carry Case should be cleaned using a soft cloth to remove debris.

15.2 - Instructions for cleaning of the Tonguetip

Before cleaning the Tonguetip ensure that it is disconnected from the Controller.

Do not expose the plug end of the Tonguetip to liquids.

- The Tonguetip should be rinsed thoroughly with water or distilled water.
- You may use the Tonguetip immediately after rinsing.
- If you do not intend to use the Tonguetip immediately after cleaning, place on a clean, dry surface and allow it to air-dry.
- In addition, the Tonguetip may be cleaned by gently rubbing the Tonguetip surfaces, including the electrodes and the plug using 70% IPA (Isopropyl alcohol) wipes.

Notes: 70% IPA wipes are available over the counter in pharmacies.

After cleaning the Tonguetip with IPA, place on a clean, dry surface to allow to air-dry for at least one minute before use.

Do not use any cleaning agent other than IPA or water.

Attempting to clean the device using a method other than that described here can damage the device and voids the warranty.

Do not scrub or wipe the surfaces forcefully.

Do not immerse the Controller, the Charger or the Tonguetip plug in water.


To avoid damaging the devices do not force dry, for example by using a hairdryer, or leaving on top of a heater.

Chapter 16 - Storage and transport

It is recommended to store the Controller and the Tonguetip in the Carry Case (model CE-401) supplied as it has been specially designed to protect them during storage and transport. It can also help maintain the cleanliness of the Tonguetip in-between uses.

It is not necessary to disconnect the Tonguetip plug from the Controller port before placing them into the Carry Case. For ease of storage, it is possible to leave the Tonguetip cable outside the case, whereby the cable enters and exits the case through the holes as indicated in this image:

The Controller and the Tonguetip should be stored in a location that is not too hot, nor too cold (between -20 and 45 degrees Celsius), and that is not excessively humid (relative humidity less than 85%).

The Controller has a shelf life of two years, meaning that it should be used for the first time within two years of it being manufactured, otherwise the battery life may deteriorate.

The environmental limits and precautions pertaining to storage and transport of the device are illustrated on the packaging labels, next to these symbols:

Storage and transport temperature limits. The lower limit is -20°C, as indicated near the lower horizontal line and the upper limit, 45°C, as indicated near the upper line.

Storage and transport humidity limits. The lower limit, 5%, as indicated near the lower horizontal line and the upper limit, 85%, as indicated near the upper line.

Keep dry.

Fragile, handle with care.

Chapter 17 - Maintenance

17.1 - Inspections

- Carry out a weekly visual inspection to ensure that the Controller has not been damaged.
- Carry out a weekly visual inspection to ensure there is no evidence of discolouration or damage to the Tonguetip.

If you have any concerns following these inspections, cease using the device and request assistance by contacting customer support, see page 27 for contact details.

17.2 - Controller service life

- The Controller has an expected service life of three years. Beyond three years, the device will still perform as intended except that the battery capacity may start to deteriorate.
- The Controller should be charged every 6 months if not in use. Not charging the battery for extended periods may reduce the capacity of the battery and shorten the service life of the device.

17.3 - Tonguetip service life

- The Tonguetip has an expected service life of three years. It is recommended not to use the Tonguetip after this time period as the electrodes may not function as effectively.
- The Tonguetip will expire after 180 cumulative hours of actual use. If an expired Tonguetip is connected to the Controller the Tonguetip Indicator will flash red and treatment sessions will not commence until a non-expired Tonguetip is connected. A new Tonguetip (Model TP-101) can be purchased from customer support. See page 27 for contact details.

Chapter 18 - Disposal

Please consider your personal safety and the environment when disposing of the Lenire system. When the device, or any of its components have reached the end of their service life, they should be disposed of in accordance with local regulations, or returned to the manufacturer (see page 27 for contact details).

No part of the Lenire system should be disposed of as domestic waste.

Section D

Troubleshooting

The following section describes problems commonly encountered while using the Lenire device, and their most common solutions. If a particular problem cannot be rectified, you should contact customer support (see page 27 for contact details).

Chapter 19 - Frequently asked questions

19.1 - What should I do if the Controller will not power on?

- Ensure that the Controller is charged (see page 11).
- Ensure the Charger is not connected. The device will not power on while the Charger is connected.
- Make sure you hold the Power Button for at least 2 seconds when attempting to turn on the device.
- It is possible that the Controller is in a state that requires a reset – connect the Charger briefly to reset the Controller and then try to turn the Controller on again once the Charger has been disconnected.
- If none of the above steps resolve the problem, please contact customer support (see page 27).

19.2 - What should I do if there is no sound coming from the Headphones?

- Ensure the Headphones are charged (see page 11).
- Ensure the Headphones are paired to the Controller. If they are not, pair the Headphones with the Controller again (see page 23).
- Check that the Headphones connect to the Controller– the Headphone indicator should illuminate green when they are connected.
- It is important to confirm that the Headphones are what the Controller is connected to (i.e. it is possible that your Controller is inadvertently connecting to another device). You can verify this by turning the Headphones off while the Controller is still on, and the Headphones indicator on the Controller should revert to pulsing blue.
- Please ensure that the Tonguetip is connected and that it has not reached expiry i.e. the Tonguetip Indicator is illuminated either green or orange, but not red.
- Ensure that the treatment session is active by pressing the play button. The Treatment Indicator will rotate around the light display when the treatment session is active.
- If there is still no sound audible through the Headphones, please contact customer support (see page 27).

19.3 - What should I do if the sound level is too low?

Try adjusting the volume level using the buttons on the side of the Controller and using the device in a quiet location. See page 17 for guidelines on adjusting the volume level.

It is likely that the sound levels will be soft, depending on the device configuration and your audiological profile.

The audio volume is calibrated to the optimum levels appropriate to your hearing thresholds during the device fitting with a trained clinician. The range of adjustment above and below this optimum level is limited in order to maximise the treatment effectiveness and minimise possible hearing damage.

19.4 - What if I find the tongue stimulation level is too low?


It is possible that there may be little or no perceived sensation from the Tonguetip (see page 16).

- Ensure that the Tonguetip Indicator illuminates green (non-pulsing) when the Tonguetip is making contact with your tongue. If it is not, please contact customer support (see page 27).
- Ensure that the Tonguetip status indicator is not flashing red. If it is then please contact the manufacturer for replacing the Tonguetip (see page 27).
- Try adjusting the stimulation level using the buttons on the side of the Controller.


19.5 - What should I do if the Controller is not charging?

If your Controller is not charging, please try the following:

- Check the mains plug on the Charger is fully 'clicked in'.

- Check if the Charger is working by connecting the Charger to the Headphones. If the red light does not illuminate, then there is possibly a fault with the Charger. Please contact customer support (see page 27).
- Retry charging the Controller - if there is still not an indication of charging, please contact customer support (see page 27).

19.6 - I tend to fall asleep during the treatment- Is this a problem?

Falling asleep during a treatment session is not advised. It is important that the Tonguetip remains in place and contacting your tongue throughout the treatment session.

19.7 - I only received one charger in the box, is that correct?

Yes. The same charger is used to charge the Controller and the Headphones.
See page 11 for details.

19.8 - What should I do if I miss a treatment session?

To maximise the effectiveness of the treatment, it is important to adhere to the treatment regimen prescribed by your clinician. If you miss a treatment session on any particular day, it is advisable that you continue treatment as normal the following day.

19.9 - My Headphones won't connect, even though the light on the Headphones indicates that it is connected

It is possible that another device, such as a mobile phone or a personal computer, is connecting to your Headphones thereby preventing the Controller from connecting to them. Try either

- Turning off any mobile phones, or other computer devices that are close by.
- Moving to a different location, away from any computer devices.

If none of these steps fixes the problem, please contact customer support (see page 27).


Chapter 20 - Bluetooth pairing / un-pairing

The Controller sends audio to the Headphones wirelessly using Bluetooth, and for this to work, the Controller and the Headphones need to be 'paired', allowing the two devices to connect automatically henceforth.


20.1 - Pairing Headphones to the Controller

The first time you try to use the Controller and Headphones together you may need to pair the devices as follows:

- Put the Headphones into pairing mode by first ensuring they are powered off, and then powering them on by briefly pressing the Headphones power button as indicated below. The light adjacent to the Headphones power button will flash blue rapidly to indicate that the Headphones are in pairing mode.
- Switch on the Controller by pressing the power button for 2 seconds.
- Place the Controller into pairing mode by pressing the play, pause and volume + buttons at the same time. All indicator lights will flash blue for 2 seconds.

- The Controller will search for and attempt to pair with the Headphones. The Headphones Indicator on the Controller will pulse rapidly during this procedure.
- If no discoverable Headphones are found within 20 minutes, the Controller will shut down. Otherwise, the Controller will connect to the Headphones and the Headphones indicator on the Controller will illuminate green.

Chapter 21 - Servicing

The device is not intended to be user-serviced. If it should fail for any reason, please return to customer support for repair (see page 27 for contact details).

This device is not intended to be modified in any way.


Chapter 22 - Device errors

If any of the following error indicators occur that prevents the use of the Lenire system, please contact customer support (see page 27 for contact details). The errors will be indicated via the lights on the Controller. All but one of the lights will turn solid red. The position of the light that is not lit will indicate the error that has occurred.


After displaying the error indicator for 5 seconds, the Controller will shut down.

Error indicators


Controller hardware failure.


Controller software error.

Internal storage error.

Tonguetip hardware failure

Tonguetip software error.

Section E

Technical information

The following section outlines technical information pertaining to the Lenire device. This information is included in order to conform to international safety standards for medical devices. There is also important information included in the following section relating to environmental operating conditions and use in the vicinity of other electrical devices such as mobile phones.

Chapter 23 - Environmental operating conditions

Lenire should only be used in environments between 5 and 40 degrees Celsius, relative humidity between 15% and 93%, and atmospheric pressure of between 700 and 1060 hPa. This precludes the use at altitudes above 3000m, unless traveling in an airplane that has the cabin pressurised.

Chapter 24 - Radio frequency (RF) specifications

The Lenire Controller contains a Bluetooth wireless radio module. Its specifications are as follows:

SPECIFICATION	DESCRIPTION
Frequency band	2.4 GHz to 2.4835 GHz
Modulation	8 DPSK, PI/4 DQPSK, GFSK
Operating range	20m to 30m
RF sensitivity	0.1% BER at -88dBm
Effective radiated power (ERP)	<4dBm (2.5 mW) Bluetooth Class 2

Chapter 25 - Guidance and declaration for electromagnetic compatibility

The Lenire system has been tested for electromagnetic performance to the following standards:

- » IEC61000-4-3 (radiated electromagnetic fields)
 - 80 MHz to 2.7 GHz at 10V/m
 - 2.7 GHz to 6 GHz at 3V/m
 - Spot frequencies per EN 60601-1-2:2015 Table 9 at up to 28 V/m
- » IEC 61000-4-2 (Electrostatic Discharge)
 - +/- 8kV Contact Discharge
 - +/- 15kV Air Discharge
- » IEC 61000-4-8 (50Hz/60Hz magnetic fields) at 30 A/m

Disturbances higher than these limits are unlikely to be encountered in the home healthcare or professional healthcare environments.

The Lenire system is intended for use in the home healthcare and professional healthcare environments only. The performance of the Lenire in environments with high electromagnetic disturbances, e.g. near magnetic resonance imaging (MRI) or high frequency surgical equipment has not been evaluated and the Lenire system should not be used in such environments.

Portable radio communications equipment (such as mobile phones) may interfere with the performance of the Lenire system (including cables as specified by Neuromod Devices Limited) and should be kept at least 30cm away from it when in use. Otherwise, degradation of the performance of this equipment could result.

The Lenire system should not be used adjacent to or stacked with other electrical equipment because it could result in improper operation.

Use of accessories and cables other than those specified for use for the Lenire system (and as provided by Neuromod Devices Limited) could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

The Lenire system has been tested for electromagnetic emissions and complies with CISPR 11 Group 1 Class B and is therefore unlikely to interfere with the operation of other electrical equipment.

Chapter 26 - Materials in direct contact with human tissue

The Lenire system does not contain phthalates, latex, human blood derivatives or cells or tissues of human or animal origin.

The electrodes of the Tonguetip contain nickel.

Note: For a complete list of materials in contact with human tissue, contact Neuromod Devices Limited.

Chapter 27 - Declaration of conformity radio equipment directive 2014/53/EU

Hereby Neuromod Devices Limited declares that the radio equipment type are Wideband Data Transmission Equipment in the 2.4 GHz ISM band (operating in the frequency range 2.4 GHz to 2.4835 GHz) per ETSI 300 328 v.2.1.1. is in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address: www.neuromoddevices.com.

Section F

Customer support contact details

General contact details

Neuromod Devices Limited,
Rainsford St,
Dublin,
D08 R2YP, Ireland.

E-mail: info@neuromoddevices.com

Web: www.neuromoddevices.com

Customer support contact details

E-mail: techsupport@neuromoddevices.com

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